




BIG FISH

Automatic Nucleic Acid Purification System

BFEX-16E

1 Pure 2 Stable 3 High-speed



Mini but powerful

Product parameter

Model	BFEX-16E
Sample processing quantity	1-16
Permanent magnetic bar (fixed)	8×2 pcs
Sample handling range	20~1000 μ l
Sample type	Animal and plant tissues, whole blood, bacteria, viruses, etc.
Overall dimensions (mm)	320×227×332 (D×W×H)
Instrument weight	12.8 kg

TEL: +86-571-56390366

Address: No. 8 Zhongqiu Road, Fuyang District,
Hangzhou City, Zhejiang Province, China

Email: support@bigfishgene.com

Website: www.bigfishgene.com