



BIG FISH

# Automatic Nucleic Acid Purification System

## BFEX-96E




1 High throughput

2 Easy

3 Stable

4 Pure



### Product parameter

Model	BFEX-96E
Sample processing quantity	96 flux/ time
Permanent magnetic bar (fixed)	96×1 pcs
Sample handling range	20~1000 $\mu$ l
Sample type	Animal and plant tissues, whole blood, bacteria, viruses, etc.
Overall dimensions (mm)	375×616×394 (D×W×H)
Instrument weight	30 kg

TEL: +86-571-56390366

Address: No. 8 Zhongqiu Road, Fuyang District, Hangzhou City, Zhejiang Province, China

Email: [support@bigfishgene.com](mailto:support@bigfishgene.com)

Website: [www.bigfishgene.com](http://www.bigfishgene.com)