



BIG FISH

Hangzhou Bigfish Bio-tech Co., Ltd.

Nucleic Acid Purification System NUETRACTION 96

BFEX-96E



Linear movement model, clear vision, 10.1-inch large color touch screen, independently designed UI interface, direct and friendly operation software.



Using wrapped deep well heating module, reduce the temperature difference between the liquid in the tube and the set temperature, improve the efficiency of lysis and elution



Fully automated and high-throughput, 1-96 samples can be processed at one run. Equipped with Bigfish extraction kit, nucleic acid extraction can be completed ultra fast



With automatic pause function of opening the door during the experiment to avoid pollution and safety problems



The instrument is equipped with air filtration and ultraviolet disinfection, which greatly reduces the risk of experimental pollution



Three kinds of intelligent magnetic absorption modes, perfectly adapted to different types of magnetic beads



Building 6, Yinhu Innovation Center, No. 9 Fuxian Road, Yinhu Street, Fuyang District, Hangzhou City, Zhejiang Province.



0571-56390366