



BIG FISH

# Real-time Fluorescent Quantitative PCR Analyzer BFQP-4850



- ✓ Patent temperature control, ultimate precision
- ✓ Smart analysis, instant access

- ✓ Modular optical path, stable signal
- ✓ Electronic hot lid, automatic pressing

## Product parameter

|                            |                                  |
|----------------------------|----------------------------------|
| Model                      | BFQP-4850                        |
| Sample volume              | 48×0.2 ml                        |
| Consumable type            | 0.2ml 8-strip tubes, single tube |
| Fluorescence channel       | 5                                |
| ROX calibration            | No need                          |
| Hot lid                    | Automatic hot lid                |
| Overall dimensions (D×W×H) | 410×300×286 mm                   |
| Instrument weight          | 13 kg                            |

TEL: +86-571-56390366

Address: No. 8 Zhongqiu Road, Fuyang District,  
Hangzhou City, Zhejiang Province, China

Email: [support@bigfishgene.com](mailto:support@bigfishgene.com)

Website: [www.bigfishgene.com](http://www.bigfishgene.com)