




BIG FISH

Thermal cycler FC-96B

- *Extremely compact*
- *Powerful temperature control*



*Small things
reveal
the big picture*

Product parameter

Model	FC-96B
Sample volume	96×0.2 ml
Monitor	Color touch screen (7 inch)
Temp. Control mode	BLOCK or TUBE
Program storage	Boundless
Gradient/Touchdown PCR	Available
Power-off protection	Available
Overall dimensions (mm)	280×187×193 (D×W×H)
Instrument weight	7 kg

TEL: +86-571-56390366

Address: No. 8 Zhongqiu Road, Fuyang District,
Hangzhou City, Zhejiang Province, China

Email: support@bigfishgene.com

Website: www.bigfishgene.com