




BIG FISH

Thermal cycler

FC-96GE

- 1 Electronic hot lid,
one-touch operation
- 2 High-end customization,
tailored to your needs
- 3 Wide temperature range,
precise temperature control



Product parameter

Model	FC-96GE
Sample volume	96×0.2 ml (Customizable)
Monitor	Color touch screen (10.1 inch)
Hot lid	Electronic hot lid
Temp. Control mode	BLOCK or TUBE
Program storage	Boundless
Gradient/Touchdown PCR	Available
Power-off protection	Available
Overall dimensions (mm)	297×353×267.5 (D×W×H)
Instrument weight	17 kg

TEL: +86-571-56390366

Address: No. 8 Zhongqiu Road, Fuyang District,
Hangzhou City, Zhejiang Province, China

Email: support@bigfishgene.com

Website: www.bigfishgene.com