








BIG FISH

Hangzhou Bigfish Bio-tech Co., Ltd.

FastCycler Thermal Cycler FC-96GE



Electronic hot cover, one key operation,
suitable for flat cover and round cover PCR tubes

-  Power-off protection: Automatically execute the remaining unfinished procedures after the power supply is restored
-  Massive storage space, program can be downloaded through USB
-  With 36°C wide gradient range, greatly convenient annealing temperature research
-  Bilingual in Chinese and English, switching freely, accurate service to domestic and foreign customers
-  Thermoelectric semiconductor technology is adopted, with accurate temperature control and rapid temperature rise and fall, up to 5 °C / s



BIG FISH
WeChat